

Analysis report No. 221027_AO3

Sample: 5-30 % CBD/ CBG tincture in HSO/OLI/SAL Oils

Sample type: oil tincture

Method: GC-FID

Sampling protocol: c.a. 130 - 250 mg sample in 5-15 ml acetone (p.a.) stirred for 20 minutes.

sample name (lot n°)	sample n°	CBD	Δ^9 -THC*	CBC	CBG	Δ^8 -THC	CBDv	THCv	CBN	CBT	CBL
		(wt.%)									
Strawberry 10% (APXE01X-HSO4-10C-ERD)	221017.28	10.64	0.15	0.18	0.08	0.04	0.03	< 0.01	< 0.01	< 0.01	< 0.01

*Regulated cannabinoid (legal limit < 1 wt.% according to ordinance 812.121.11 of Swiss Federal Department of Home Affairs)

N.B.: Indicated values account for total amount of cannabinoids after complete decarboxylation of the material.



Bereket Tesfazghi

Chemist



Dr Jebriil HADI

Head of Laboratory