

Analysis report No. 221027_AO3

Sample: 5-30 % CBD/ CBG tincture in HSO/OLI/SAL Oils

Sample type: oil tincture

Method: GC-FID

Sampling protocol: c.a. 130 - 250 mg sample in 5-15 ml acetone (p.a.) stirred for 20 minutes.

| sample name (lot n°) | sample n° | CBD | Δ^9 -THC* | CBC | CBG | Δ^8 -THC | CBDv | THCv | CBN | CBT | CBL |
|--|-----------|--------|------------------|------|------|-----------------|------|--------|--------|--------|--------|
| | | (wt.%) | | | | | | | | | |
| Strawberry 20% (APXE01X-HSO4-20A-ERD) | 221017.31 | 20.74 | 0.16 | 0.19 | 0.08 | 0.05 | 0.06 | < 0.01 | < 0.01 | < 0.01 | < 0.01 |

*Regulated cannabinoid (legal limit < 1 wt.% according to ordinance 812.121.11 of Swiss Federal Department of Home Affairs)

N.B.: Indicated values account for total amount of cannabinoids after complete decarboxylation of the material.



Bereket Tesfazghi

Chemist



Dr Jebriil HADI

Head of Laboratory