

Analysis report No. 221027_AO4

Sample: 5-30 % CBD/ CBG tincture in HSO/OLI/SAL Oils

Sample type: oil tincture

Method: GC-FID

Sampling protocol: c.a. 130 - 250 mg sample in 5-15 ml acetone (p.a.) stirred for 20 minutes.

sample name (lot n°)	sample n°	CBD	Δ ⁹ -THC	* свс	СВG	Δ ⁸ -THC	CBDv	THCv	CBN	СВТ	CBL
	(wt.%)										
Natural 30% (APXE01X-OLI4-30A)	221024.14	30.92	0.18	0.21	0.24	0.05	0.10	< 0.01	< 0.01	< 0.01	< 0.01

N.B.: Indicated values account for total amount of cannabinoids after complete decarboxylation of the material.

Bereket Tesfazghi

Chemist

Dr Jebril HADI

Head of Laboratory

^{*}Regulated cannabinoid (legal limit < 1 wt.% according to ordinance 812.121.11 of Swiss Federal Department of Home Affairs)